

STOOL COLLECTION INSTRUCTIONS FOR OVA AND PARASITE EXAMINATION


HOW TO COLLECT A STOOL SAMPLE?

Your doctor has ordered a laboratory test that requires stool sample collection.

Please read the instructions for proper stool sample collection.

From the laboratory you will need:


- White transport container
- Brown-top container with spoon and special preservative
- Stool collection paper



- 1 Label the container in advance.




- 2 Open the stool collection paper following the arrows indicated on the illustration.



- 3 Place the stool collection paper at the back of the toilet bowl as shown in the illustration. Avoid wetting the paper.

- 4 Remove the cap from the brown-top container. Using the attached spoon, collect stool from 2-3 different areas and place the material into the same container. Tightly close the cap.




- 5 Collect the specimen daily for 3 consecutive days in the same container.

- 6 After specimen collection, remove the paper from both sides of the toilet seat and dispose of it in the toilet.

Allow a few seconds for the paper to absorb water, then flush.

- 7 For delivery to the laboratory, place the brown-top container with collected specimen into the white transport container.



Note: The specimen is stable for an extended period of time and does not require special storage.